

## STR Authentication Report For Cell Line

1. Cell Name: HGC-27-Y01-ZKY-ZZK

2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 21- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer.

3. Sample gene genotype analysis result:

Genotype analysis results of STR and Amelogenin loci in cells						
STR Loci	STR profile for Sample cell			STR profile for Cell Bank cell		
	Cell name: HGC-27-Y01-ZKY-ZZK			Cell name: HGC-27		
	Allele1	Allele2	Allele3	Allele1	Allele2	Allele3
D5S818	12	12		12	12	
D13S317	10	11		10	11	
D7S820	11	12	13.0	11	12	13.0
D16S539	10	11		10	11	
VWA	14	14		14	14	
TH01	9	9		9	9	
AMEL	X	X		X	X	
TPOX	8	8		8	8	
CSF1PO	12	12		12	12	

Note: The above figure shows the comparison results of nine genotype loci for detecting cells and cells from the cell bank.

4. Conclusion: The STR results of **HGC-27** cells are consistent with the genotypes of **HGC-27** cell lines in **DSMZ** database, the cell ID corresponded to **RCB0500** , and the STR results **completely matched**. In the test, **Multi allele was found** in this cell line.

## 5. Attached Image

