

STR Authentication Report For Cell Line

1. Cell Name: Jurkat-CloneE6-1

2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 21- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer.

3. Sample gene genotype analysis result:

Genotype analysis results of STR and Amelogenin loci in cells						
STR Loci	STR profile for Sample cell			STR profile for Cell Bank cell		
	Cell name: Jurkat-CloneE6-1			Cell name: Jurkat, Clone E6-1		
	Allele1	Allele2	Allele3	Allele1	Allele2	Allele3
D5S818	9	9		9	9	
D13S317	8	12		8	12	
D7S820	8	12		8	12	
D16S539	11	11		11	11	
VWA	18	18		18	18	
TH01	6	9.3		6	9.3	
AMEL	X	Y		X	Y	
TPOX	8	10		8	10	
CSF1PO	11	12		11	12	

Note: The above figure shows the comparison results of nine genotype loci for detecting cells and cells from the cell bank.

4. Conclusion: The STR results of **Jurkat, Clone E6-1** cells are consistent with the genotypes of **Jurkat, Clone E6-1** cell lines in **DSMZ** database, the cell ID corresponded to **TIB-152** , and the STR results **completely matched**. In the test, **Multi allele was found** in this cell line.

5. Attached Image

