

STR Authentication Report For Cell Line

1. Cell Name: KYSE30
2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 21- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer.
3. Sample gene genotype analysis result:

Genotype analysis results of STR and Amelogenin loci in cells						
STR Loci	STR profile for Sample cell			STR profile for Cell Bank cell		
	Cell name: KYSE30			Cell name: KYSE-30		
	Allele1	Allele2	Allele3	Allele1	Allele2	Allele3
D5S818	11	11		11	11	
D13S317	9	9		9	9	
D7S820	11	12		11	12	
D16S539	10	12		10	12	
VWA	16	19		16	18	19.0
TH01	9	9		9	9	
AMEL	X	X		X	X	
TPOX	9	9		9	9	
CSF1PO	10	10		10	10	

Note: The above figure shows the comparison results of nine genotype loci for detecting cells and cells from the cell bank.

4. Conclusion: The STR results of **KYSE-30** cells are consistent with the genotypes of **KYSE-30** cell lines in **DSMZ** database, the cell ID corresponded to **351** , and the STR results **basically matched**. In the test, **No multi allele was found** in this cell line.

5. Attached Image

