

STR Authentication Report For Cell Line

1. Cell Name: MHCC97-H
2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 21- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer.
3. Sample gene genotype analysis result:

Genotype analysis results of STR and Amelogenin loci in cells						
STR Loci	STR profile for Sample cell			STR profile for Cell Bank cell		
	Cell name: MHCC97-H			Cell name: MHCC97-H		
	Allele1	Allele2	Allele3	Allele1	Allele2	Allele3
D5S818	12	13		12	13	
D13S317	8	8		8	11.3	
D7S820	10	10		10	10	
D16S539	12	12		12	12	
VWA	14	14		14	14	
TH01	9	9		9	9	
AMEL	X	Y		X	Y	
TPOX	8	8		8	8	
CSF1PO	11	13		11	13	

Note: The above figure shows the comparison results of nine genotype loci for detecting cells and cells from the cell bank.

4. Conclusion: The STR results of **MHCC97-H** cells are consistent with the genotypes of **MHCC97-H** cell lines in **EXPASY** database, the cell ID corresponded to **CVCL_4972** , and the STR results **basically matched**. In the test, **Multi allele was found** in this cell line.

5. Attached Image

