

## STR Authentication Report For Cell Line

1. Cell Name: SK-HEP-1
2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 21- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer.
3. Sample gene genotype analysis result:

Genotype analysis results of STR and Amelogenin loci in cells						
STR Loci	STR profile for Sample cell			STR profile for Cell Bank cell		
	Cell name: SK-HEP-1			Cell name: SK-HEP-1		
	Allele1	Allele2	Allele3	Allele1	Allele2	Allele3
D5S818	10	13		10	13	
D13S317	8	12		8	12	
D7S820	8	11		8	11	
D16S539	12	12		12	12	
VWA	14	17		14	17	
TH01	7	9		7	9	
AMEL	X	X		X	X	
TPOX	9	9		9	9	
CSF1PO	11	12		11	12	

Note: The above figure shows the comparison results of nine genotype loci for detecting cells and cells from the cell bank.

4. Conclusion: The STR results of **SK-HEP-1** cells are consistent with the genotypes of **SK-HEP-1** cell lines in **DSMZ** database, the cell ID corresponded to **141** , and the STR results **completely matched**. In the test, **No multi allele was found** in this cell line.

### 5. Attached Image

