

STR Authentication Report For Cell Line

1. Cell Name: U-2 OS
2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 21- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer.
3. Sample gene genotype analysis result:

Genotype analysis results of STR and Amelogenin loci in cells						
STR Loci	STR profile for Sample cell			STR profile for Cell Bank cell		
	Cell name: U-2 OS			Cell name: U-2 OS		
	Allele1	Allele2	Allele3	Allele1	Allele2	Allele3
D5S818	11	11		11	11	
D13S317	13	13		13	13	
D7S820	11	12		11	12	
D16S539	11	12		11	12	
VWA	14	18		14	18	
TH01	6	9.3		6	9.3	
AMEL	X	X		X	X	
TPOX	11	12		11	12	
CSF1PO	13	13		13	13	

Note: The above figure shows the comparison results of nine genotype loci for detecting cells and cells from the cell bank.

4. Conclusion: The STR results of **U-2 OS** cells are consistent with the genotypes of **U-2 OS** cell lines in **DSMZ** database, the cell ID corresponded to **HTB-96**, and the STR results **completely matched**. In the test, **No multi allele was found** in this cell line.

5. Attached Image

